3D time lapse imaging and quantitative analysis of the active migration of human vascular endothelial cells into a multilayered cell sheet with Yokogawa's confocal image cytometer CQ1

Toshiaki Endo¹; Yoichi Nomura¹, Michael Schell², Romina Dintner², Aki Ota¹ ¹Yokogawa electric corporation, 2-3 Hokuyodai Kanazawa, Japan ²Cenibra GmbH, Große Str 17 D-49565 Bramsche, Germany

Introduction

Time lapse confocal imaging has been an essential method to investigate the 3D dynamic behaviors of cells in tissue cultures. For long-term live cell imaging, it is critical to reduce phototoxic damage to the cells caused by repeated laser scanning. Yokogawa CSU (confocal scanner unit) is a confocal unit using a microlens-enhanced dual Nipkow disk confocal optical system, which has been shown to be less harmful to living cells compared to conventional single beam scanning devices. The CQ1 is an all-in-one confocal quantitative imaging cytometer based on the CSU. Here we report the 3D time lapse live cell imaging in a multilayered cell sheet using CQ1.

Results

1. Dynamic migration and network formation of GFP-HUVECs captured by 3D time lapse imaging



Methods

Five-layered myoblast cell sheets were constructed from human skeletal muscle myoblasts (HSMM) and human skeletal muscle fibroblasts (HSMF) by using a thermo-responsive surface and stamp system ⁽¹⁾⁽²⁾. HSMMs and HSMFs were labeled with CellTracker[™] Orange. Human umbilical vein endothelial cells (HUVEC) expressing GFP (GFP-HUVEC) were overlaid by the cell sheet and co-cultured. Time lapse imaging (67 hours, 30 min interval, 40x objective lens, 49 fields) was performed by CQ1 equipped with an internal incubation chamber to regulate culture environment.





Fig. 1-1. 3D images of the cell sheet at the beginning of and after 51 hours incubation. Images were reconstructed from Z stack images of the field indicated by the yellow frame in the large field stitched image in Method fig.2.









Method fig. 1. The fivelayered cell sheet

Method fig. 2. Large field stitched 3D image acquired by CQ1 7 x 7 fields, Green; HUVEC (GFP), Red; Cell Sheet (CellTrackerTM Orange)

YOKOGAWA Core Competence

Microlens enhanced dual Nipkow disk scanning

 \rightarrow high-speed, low photo toxicity and low photo bleaching

A Nipkow spinning disk and a second spinning disk with microlens rapidly scan the field of view with about 1,000





51hr 0hr 34hr 17hr

Fig. 1-2. Migration of the GFP-HUVECs into the cell sheet.

Single slice images showing the migration of HUVECs into upper layers. (*Rows, from top to bottom*) Single slice images of layers 3, 2, 1 and corresponding Y-Z plane images of the cell sheet. (Columns, from left to right) Images acquired at 0, 17, 34 and 51 hr incubation. The image filed is the same as fig. 1-1.

2. Quantification of the migration of GFP-HUVECs into the five-layered cell sheet



Summary & Conclusions

GFP-HUVECs dynamically migrated upward into the five-layered cell sheet constructed from HSMMs and HSMFs.

Confocal Quantitative Image Cytometer CQ1

üEasy to use to measure the feature data of cells üCell clusters are directly measured by 3D imaging

üGently 3D image acquisition

üBench-top and no need for darkroom



- The GFP-HUVECs formed a reticulate network in the horizontal plane in the middle layers.
- Subscription Long-term 3D time lapse imaging by CQ1 revealed a dynamic process of the active migration and the formation of the cellular network in the multilayered cell sheet.
- **Ø** CQ1 would be a powerful research tool in tissue engineering as well as regenerative medicine and drug screening.

Acknowledgements

This study was conducted under the supervision of Dr. Nagamori, Osaka University and with Yokogawa Electric Corporation's responsibility.

Reference

(1) Kino-oka et al., J. Biosci. Bioeng., 113, 128-131 (2012) (2) Nagamori et al., Biomaterials, 34, 662-668 (2013)