

RNAGEM V

Effortless, simultaneous, nucleic acid extraction from DNA and RNA viruses ... in minutes not hours.

RNAGEM V utilizes MicroGEM's novel thermophilic proteinase and buffer system to extract nucleic acids from viruses in less than 15 minutes without the need for purification.

Sample types: nasal-pharyngeal swabs in saline solution, nasal swabs, buccal swabs & sputum/saliva

Downstream applications: PCR, RT-PCR, RT-qPCR, LAMP, RT-LAMP

Key advantages:

- Nucleic acid extraction in less than 15 minutes
- Simultaneous extraction from DNA and RNA viruses
- Single-tube extraction
- No ionic detergents or chaotropic salts
- High nucleic acid recovery – Minimal loss of nucleic acids during extraction
- No magnetic beads – No spin-columns
- Flexibility - suitable for low-throughput to high-throughput extraction with a single protocol
- Easily automated using standard liquid handling solutions
- Minimal plasticware required - Reduced waste and supply chain issues

RNAGEM V is for research use only (RUO), not intended for in-vitro diagnostic (IVD) use.

Image shows 50 and 1,000 reaction kits.



Typical Workflow

RNAGEM V protocols can be scaled up or down as needed.



1. Prepare sample.

2. Mix sample and *RNAGEM V* reagents.



3. Heat mixture at 75° C to activate *RNAGEM* enzyme.

4. Denature *RNAGEM* enzyme at 95° C.



5. This solution now contains viral nucleic acid ready for PCR based applications.

