

User Guide

High Sensitivity (N1) Cartridge Kit (C105105/C105205)

A. Specifications

| Specification | Description |
|------------------|-------------|
| DNA Sizing Range | 20-5,000 bp |
| L.O.D | 5 pg/μL |
| Resolution* | 4-10 bp |
| Sample Number | 100 runs |
| Shelf Life | 4 months |

*Best resolution is determined by the 15-622 DNA Size Marker (C109200).

B. Kit Components and Storage Conditions

| Item | Storage Condition |
|---|--|
| High Sensitivity Cartridge (C105105/C105205) | 4-27°C (Do Not Freeze) |
| 20-1,000 bp Alignment Marker (C109100-500A, 100 μL) | Short-Term (≤ 6 months): 4-27°C Long-Term (> 6 months): -20°C |
| Separation Buffer (C104406, 50mL) | 4-27°C |
| Dilution Buffer (C104405, 15mL) | 4-27°C |
| Mineral Oil (C104404, 8mL) | 4-27°C |

C. Cartridge Unpacking Preparation

New cartridge must pass HV check and calibration before use. Please follow unpacking guide to unpack and use C109100-500A Alignment to do calibration.

D. Buffer and Marker Preparation

Qsep₁₀₀

Buffers: add to the groove

Qsep₁ 8-well tray

Qsep₁ 12-well tray

Qsep₁-Plus 16+3 well tray

Separation Buffer S: 400 μl Separation Buffer
Park P: 200 μl d₂H₂O
Clean/Wash Buffer C: 400 μl d₂H₂O

E. Sample Preparation

Sample Volume Requirements

0.2 mL Tube: 20 μL Micro Vial (C104250): 2 μL
0.1 mL Tube: 10 μL 16+3-Well Tube (C104254): 10 μL

Recommended Sample Concentration

Fragment sample: 5-200 pg/μl

***NOTE: When fragment sample concentration is over 500 pg/μl, dilute sample 10X with 0.5X dilution buffer.**

Smear sample (NGS library/cfDNA): 50-1000 pg/μl

***NOTE: When smear sample concentration is over 1000 pg/μl, dilute sample 10X with 0.5X dilution buffer.**

***NOTE: If sample is eluted in water, add dilution buffer to make the sample into 0.5X dilution buffer condition.**

E-1. Sample size within range from 20 bp to 1000 bp

Markers required:

20 bp-1000 bp Alignment Marker (C109100): 20 μl
20X diluted 15-622bp Size Marker (C109200): 20 μl

| Materials | Volume (μl) |
|----------------------------------|-------------|
| Size Marker (15-622bp) (C109200) | 1 |
| Dilution Buffer (C104402) | 19 |
| Total | 20 |

E-2. Sample size within range from 20 bp to 5000 bp

Markers required:

20 bp-5000 bp Alignment Marker (C109102): 20 μl
10X diluted 50-3000bp Size Marker (C109301): 20 μl

| Materials | Volume (μl) |
|-----------------------------------|-------------|
| Size Marker (50-3000bp) (C109301) | 2 |
| Dilution Buffer (C104402) | 18 |
| Total | 20 |

***NOTE: If sample is eluted in nuclease-free water, using dilution buffer to make the sample into 0.5X dilution buffer condition is suggested.**


Contact Information:

Company Name: Bioptic Inc.
Address: (23141) 5F., No.6, Ln. 130, Minquan Rd., Xindian District, New Taipei City, Taiwan (R.O.C)
Tel: +886-2-2218-8726, Fax: +886-2-2218-8727, E-mail: service@bioptic.com.tw


User Guide

High Sensitivity (N1) Cartridge Kit (C105105/C105205)


F. Software Operation



1. Place the sample and select corresponding position, and then input sample information.



2. Select the alignment marker and the analytic method "cfDNA" in Method Selector.



***NOTE: Based on sample concentration to adjust injection condition**

| Sample concentration | High (4kV, 20s) | Regular (6kV, 20s) | Low (8kV, 20s) |
|----------------------|-----------------|--------------------|----------------|
| Fragment DNA | >200 pg/μl | 10~200 pg/μl | <10 pg/μl |
| Smear DNA | > 1000 pg/μl | 50-1000 pg/μl | <50 pg/μl |

3. Sample Duration: adjust the sample injection time to increase/decrease injection amount.

***NOTE: Do not set the injection time over 20 sec.**


4. Runs: reputation time.

5. Separation Duration: adjust the duration to extend/reduce the separation time.

***NOTE: Step 3-5 are optional.**

6. Input the result name for result file and select the tab below to bring information into result name.

7. Click "Para" . Choose to use (a) reference or (b) create size marker to do the calculation.



***NOTE: When using function "Create Size Marker", select the size marker you use. e.g. 20-1k is paired with C109200, and 20-5k is paired with C109301.**

8. Click "Run" to start analysis.

Contact Information:

Company Name: BiOptic Inc.

Address: (23141) 5F., No.6, Ln. 130, Minquan Rd., Xindian District, New Taipei City, Taiwan (R.O.C)



Tel: +886-2-2218-8726, Fax: +886-2-2218-8727, E-mail: service@bioptic.com.tw

User Guide


High Sensitivity (N1) Cartridge Kit (C105105/C105205)

G. Result and Application


• Alignment Marker & Size Marker



• Cell-Free DNA



• Cell-Free DNA with genomic DNA



H. Troubleshooting

Please ensure that the system is working well, and the operation follows the instructions first.

Sometimes, there will be some residues left in cell-free DNA sample after extraction. These residues might cause unstable current at sample injection or separation steps. Here is a list of solutions to help fix the occurrence.

1. Use dilution buffer or nuclease-free water to dilute the sample.
***NOTE: Final buffer concentration in 0.5X Dilution**
2. Centrifuge the sample for a while to make the residues accumulate at the bottom of the tube.
3. Insert a "T-purge-120" method between several sample runs.e.g. insert 1 run "T-Purge-120" every 5-10 sample runs.

| Sequence | | | | | | | | Open | Save | Save As |
|----------|-----------------|-------------|-----------------|------|---------------------|-------------|-----------|------|--------|---------|
| SN | Sample Position | Method | Sample Duration | Runs | Separation Duration | Result Name | Para | Add | | |
| 1 | A-01,A-0 | cdDNA | 20 | 1 | 200 | Test00 | Sample ID | | Insert | |
| 2 | | T-Purge-120 | 0 | 1 | 0 | Sample ID | | | Delete | |
| 3 | A-06,A-0 | cdDNA | 20 | 1 | 200 | Test01 | Sample ID | | Up | |
| | | | | | | | | | Down | |

e.g. Insert 1 run "T-Purge-120" every 5-10 sample runs

Micro Vial

00:47

I. Cartridge Discard

Please wear gloves before discarding cartridge.

Gel reservoir



Cartridge tip

1. Bend the cartridge tip.
2. Open the cap on gel reservoir and remove the inner cap.
3. Pour the gel into the chemical waste container.
4. Cartridge can be thrown into the bin.

Contact Information:

Company Name: BiOptic Inc.
Address: (23141) 5F., No.6, Ln. 130, Minquan Rd., Xindian District, New Taipei City, Taiwan (R.O.C)
Tel: +886-2-2218-8726, Fax: +886-2-2218-8727, E-mail: service@bioptic.com.tw