



Yeastern Biotech Co., Ltd



# T&A™ Cloning Vector Kit Protocol

## Product Use Limitation & Warranty

This product is intended to be used for life science research only. It has not been approved for drug or diagnostic purpose. YEASTERN's products should not be resold, modified for resale, or used to manufacture commercial products without written approval by YEASTERN. YEASTERN guarantees the performance of all products in the manner described in our protocol. The purchaser must determine the suitability of the product for its particular use. Should any product fail to perform satisfactorily due to any reason other than misuse, YEASTERN will replace it free of charge.

Ver. L0514

No part of these protocols may be reproduced in any form or by any mean, transmitted, or translated into a machine language without the permission of YEASTERN BIOTECH CO., LTD.

Address: 6F-3, 23 Lane 169, Kang Ning St., Shijr, Taipei, 22180 Taiwan.  
Tel: +886-2-2695-3922 Fax: +886-2-2695-3979  
Email: [yeastern@yeastern.com.tw](mailto:yeastern@yeastern.com.tw) Website: [www.yeastern.com](http://www.yeastern.com)

Copyright© 2012 All rights reserved. Yeastern Biotech Co., Ltd.


Copyright© 2012 All rights reserved. Yeastern Biotech Co., Ltd.

Cat. No.  
FYC001-20P  
FYC002-20P

# Manual Index

1. Vector map of T&A™ Cloning Vector	2
2. DNA sequence of multiple cloning sites in T&A™ Cloning Vector	3
3. Protocol for ligation using T&A™ Cloning Vector	3
4. Protocol for colony PCR	4
5. Product components	5
6. Restriction enzyme sites of T&A™ Cloning Vector	6
7. Suggestions and notes	8
8. Troubleshooting	9

## 1. Vector map of T&A™ Cloning Vector



**Figure 1. Map and sequence reference points of T&A™ Cloning Vector**

\* Before the insert is incorporated into T&A™ Cloning Vector, there is only one *Hind*III site and no *Bgl*II site present in the multiple cloning region. After the incorporation, the 3'-A overhangs at both end of the insert will complement with the 3' overhanging T at the terminals of the vector and generate one additional *Hind*III site and one new *Bgl*II site. This merit of T&A™ vector makes cloning more economical and convenient.

Multiple cloning region	434 to 490
<i>LacZ</i> start codon	511
<i>LacZ</i> operator	531 to 548
<i>LacZ</i> gene	511 to 149
<i>Amp<sup>r</sup></i> gene	2528 to 1671
T7 promoter	402 to 439
M13 forward primer	359 to 375
M13 reverse primer	507 to 528
$\beta$ -lactamase coding region	1524 to 2528
<i>Lac</i> operon sequences	239 to 398, 488 to 725

## 2. DNA sequence of multiple cloning region in T&A™ Cloning Vector

301	TACGCCAGCT ATGCGGTCGA	GGCGAAAGGG CCGCTTTC	GGATGTGCTG CCTACACGAC	CAAGGCGATT GTTCCGCTAA	AAGTTGGGTA TTCACCCAT
	<b>M13 Forward Primer</b>				
351	ACGCCAGGGT TGCGGTCCCA	TTTCCAGTC AAAGGGTCAG	ACGACGTTGT TGCTGCAACA	AAAACGACGG TTTTGCTGCC	CCAGTGAATT GGTCACTTAA
	<b>T7 Promoter</b>				
401	GTAATACGAC CATTATGCTG	TCACTATAGG AGTGATATCC	GCGAGCTCGG CGCTCGAGCC	TACCCGGGCG ATGGGCCCGC	AATTTCCAAGC TTAAGGTTCG
		<i>Bgl</i> II	<i>Bam</i> HI	<i>Xba</i> I	<i>Sal</i> I
451	TT AA	AGATCTGGAT TCTAGACCTA	CCCCCTCTAGA GGGGAGATCT	GTCGACCTGC CAGCTGGACG	AGGCATGCAA TCCGTACGTT
		<i>Hind</i> III		<i>Pst</i> I	
493	CGTTGGCGTA GCAACCGCAT	ATCATGGTCA TAGTACCAGT	TAGCTGTTTC ATCGACAAG	CTGTGTGAAA GACACACTTT	TTGTTATCCG AACATAGGC
	<b>M13 Reverse Primer</b>				

**Figure 2. Multiple cloning site sequences of T&A™ Cloning Vector**

## 3. Protocol for ligation using T&A™ Cloning Vector

1. Centrifuge T&A™ Cloning Vector and PCR DNA samples briefly to collect contents at the bottom of the tubes.
2. Vortex the ligation buffers vigorously before use.
3. Set up the following items as described below :

	Standard Control	Positive Control
Ligation Buffer A	1 $\mu$ l	1 $\mu$ l
Ligation Buffer B	1 $\mu$ l	1 $\mu$ l
T&A™ cloning vector	2 $\mu$ l	2 $\mu$ l
PCR product	X $\mu$ l	****
yT <sub>4</sub> DNA ligase	1 $\mu$ l	1 $\mu$ l
Control DNA	****	3 $\mu$ l
Add deionized water to a final volume of 10 $\mu$ l		


4. Mix the reactions by pipetting.
5. Incubate the reactions for 5 to 15 min at 22°C. Alternatively, if higher transformation efficiency is needed, incubate the reactions overnight at 4°C.
6. Transformation with appropriate competent cells according to user's downstream applications.

## 4. Protocol for colony PCR

1. Prepare 25 µl of PCR reaction buffer in a microfuge tube as described below. Pick a single colony with a sterile toothpick and directly swirl it into PCR reaction buffer as the template.

PCR premixed buffer (O <sup>3</sup> in1 DNA polymerase premix, FYT201-100P)	23 µl
M13-F (10 µM)	1 µl
M13-R (10 µM)	1 µl


2. Set up PCR program



<sup>+++</sup>The time duration at 72°C is set according to the length of insert DNA.  
(In general, DNA polymerase can synthesize 1 Kb of DNA in 1 minute.)

3. Run on 1% agarose gel

Example : Using the control DNA provided in the T&A™ cloning vector kit as the insert DNA. The colony PCR result is shown as below.



\*According to different primer designs, the size of colony PCR product seen on the agarose gel will change.  
In this colony PCR, the size of the PCR product on the gel is larger than the insert DNA size by about 150 bp.

## 5. Product components

T&A™ Cloning Kit	Cat. No. FYC001-20P	
Components	Concentration	Volume
T&A™ Cloning Vector (20 reactions)	25 ng/µl	40 µl
Control Insert DNA	10 ng/µl	10 µl
yT <sub>4</sub> DNA Ligase	2 U/µl	20 µl
10x Ligation Buffer A	----	50 µl
10x Ligation Buffer B	----	50 µl
Forward Primer (M13-F)	10 µM	50 µl
Reverse Primer (M13-R)	10 µM	50 µl
Storage Condition: -20°C		

T&A™ Cloning Vector	Cat. No. FYC002-20P	
Components	Concentration	Volume
T&A™ Cloning Vector (20 reactions)	25 ng/µl	40 µl
Control Insert DNA	10 ng/µl	10 µl
Storage Condition: -20°C		

## Suggestions

1. **Avoid multiple freeze-thaw cycles** and exposure to frequent temperature changes by preparing single-use aliquots of the ligation buffers.
2. *Pfu* DNA polymerase possesses proofreading activity; it does not have the terminal transferase-like activity like Taq DNA polymerase. Ligation reactions using *Pfu* amplified DNA containing no A-tails will result in no positive colonies.
3. Methods for increasing the ligation efficiency :
  - A. A-tailing:
 

purified PCR product	X µl
10X PCR buffer	10 µl
10 mM dATP	2 µl
Taq	1 µl

- i. Add deionized water to a final volume of 100  $\mu$ l.
  - ii. Incubate at 72°C for 1 hr.
- Purify the A-tailed DNA and use in the ligation reaction.

- B.** If the maximum of transformants is required, incubate the reaction overnight at 4°C.
- C.** The optimal efficiency can be achieved by using a 1:3 molar ratio of vector DNA to the insert DNA.
- D.** Use competent cells with higher efficiency such as ECOS™ (>10<sup>8</sup> cfu/ $\mu$ g DNA) series for transformation.

## 6. Restriction enzyme sites of T&A™ Cloning Vector

Enzyme	Position	Enzyme	Position	Enzyme	Position	Enzyme	Position	Enzyme	Position
Aafl	2664	AspEI	1742	Cfr10I	1822	MamI	457	SspI	2546
Acc65I	430	AvaI	434	Drall	2718	NarI	237	XbaI	466
AccI	473	BanII	428	Eam1105I	1742	NdeI	185	XmaI	434
AcsI	441	BamHI	458	Ecl136II	426	PstI	482	XmnI	2341
AflIII	849	Bcgl	2281	Eco109I	2718	SacI	428		
AhdI	1742	Bpml	1812	EcoRI	441	SalI	472		
AlwNI	1265	BsaBI	457	HincII	474	SapI	733		
ApoI	441	BsaI	1803	HindII	474	Scal	2222		
Asp700	2341	BspMI	485	KasI	236	Smal	436		
Asp718	430	BsrFI	1822	KpnI	434	SphI	488		

## Restriction enzymes that DO NOT cut T&A™ Cloning Vector

Aafl	BbsI	BsNI	BstBI	Eco47III	MroI	NsiI	Ppu10I	SpeI
AccIII	BclI	BsmFI	BstEII	EcoNI	MscI	NspV	PpuMI	SspBI
AflII	BfrI	BsmI	Bsu36I	EcoRV	MunI	Pacl	RsrII	StyI
AgeI	BlnI	Bsp120I	CelII	EspI	NaeI	PaeR7I	SacII	Swal
Apal	Bpu1102I	BspDI	Clal	HpaI	NcoI	PfMI	SexAI	Tth111I
AspI	BpuAI	BspEI	Csp45I	KspI	NgoMI	PinAI	SfiI	Van91I
AsuII	BsaAI	BsrGI	Drall	MfeI	NheI	PmaCI	SfuI	XcmI
AvrII	BseAI	BssHII	EagI	MluI	NotI	PmeI	SgrAI	XhoI
BbrPI	BsgI	Bst1107I	EcXI	MluNI	NruI	PmlI	SnaBI	XmaIII

## Restriction enzymes that cut T&A™ Cloning Vector more than 2 times

Name	Position	Name	Position	Name	Position
AcI	34 sites	CfoI	17 sites	MboI	16 sites
AluI	17 sites	Csp6I	169, 431, 2221	MbolI	7 sites
Alw44I	178, 1163, 2409	DdeI	6 sites	McrI	5 sites
AlwI	10 sites	DpnI	16 sites	MnlI	13 sites
AosI	259, 1964	DpnII	16 sites	MseI	13 sites
ApaLI	178, 1163, 2409	DraI	1608, 1627, 2319	MslI	1994, 2153, 2512
AseI	620, 679, 1914	DrdI	98, 957	MspA1I	6 sites
AsnI	620, 679, 1914	DsaV	12 sites	MspI	13 sites
AspHI	5 sites	EaeI	389, 688, 2130	MvaI	5 sites
Avall	1880, 2102	EarI	297, 733, 2537	Mvnl	10 sites
Avll	259, 1964	Eco57I	1397, 2409	MwoI	13 sites
		EcoRII	5 sites	NciI	7 sites
BanI	4 sites	Esp3I	4, 46	NdeII	16 sites
BbvI	12 sites	Fnu4HI	19 sites	NlaIII	11 sites
Bfal	4 sites	FnuDII	10 sites	NlaIV	11 sites
BglI	252, 1862	FokI	5 sites	NspI	42, 488, 853
BmyI	5 sites	FspI	259, 1964	PleI	5 sites
BsaHI	237, 2279, 2661	HaeII	240, 727, 1097	Psp1406I	1968, 2341
BsaJI	4 sites	HaeIII	11 sites	PvuI	280, 2112
BsaWI	1055, 1202, 2033	HgaI	4 sites	PvuII	309, 673
BsEI	5 sites	HgiAI	5 sites	RcaI	1569, 2577, 2682
BsHKAI	5 sites	HhaI	17 sites	RsaI	170, 432, 2222
BsYI	6 sites	HindIII	(448), 490	Sau3AI	16 sites
BsI	6 sites	Hinfl	7 sites	Sau96I	6 sites
BsmAI	4 sites	HinPI	17 sites	ScrFI	12 sites
Bsp1286I	5 sites	HpaII	13 sites	SfaNI	8 sites
BspHI	1569, 2577, 2682	HphI	7 sites	Sfcl	5 sites
BspWI	13 sites	HPY188I	10 sites	Snol	178, 1163, 2409
BsrBI	541, 782, 2583	Ital	19 sites	TaqI	473, 949, 2393
BsrDI	1803, 1977	Ksp632I	297, 733, 2537	TfiI	684, 824
BsrI	11 sites	MaeI	4 sites	Thal	10 sites
BstNI	5 sites	MaeII	5 sites	Tru9I	13 sites
BstUI	10 sites	MaeIII	11 sites	Tsp509I	8 sites
BstYI	7 sites			XhoII	7 sites

## 7. Suggestions and notes

1. For blue-white selection, plating should be done onto agar plates containing 20-50 µg/ml ampicillin plus 0.1-0.5 mM IPTG and 40-60 µg/ml X-gal. Alternatively, the IPTG and X-gal can be spread evenly onto a 20-50 µg/ml ampicillin plate at least 2 hr prior to plating. After the colonies grow for at least 14 hr, the plates can be put at 4°C and the blue color will continue to develop.
2. When coupled with the colony PCR technique, clones can be screened easily and precisely.
3. Multiple freeze-thawing DOES NOT affect the quality of T&A™ Cloning Vector, but exposure to frequent temperature changes will degrade ATP resulting in poor ligation.
4. For questions not addressed here, please visit our web site for details. <http://www.yeastern.com>

## 8. Troubleshooting

Symptom	Comments	Suggestions
No colonies	Problems in transformation or use low-efficiency competent cells.	Self-ligated vector should only yield < 50 colonies, while with A-tailed insert DNA control the transformation should yield > 500 colonies. If there is no colony on the plate, use high-efficiency competent cells to overcome the problem.
Less than 100 transformed colonies grow on selection plate. When checking the colonies with colony PCR using M13F&R primers, false inserts of 0.15 and 1.8 kb fragments are obtained.	<ol style="list-style-type: none"> <li>1. Improper molar ratio of the vector DNA to the insert DNA.</li> <li>2. Bad A-tailed insert DNA.</li> <li>3. The 0.15kb fragment is derived from reaction of self-ligation of the vector itself. The 1.8 kb fragment is derived from the residue in production process.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check the A-tailed insert DNA quality.</li> <li>2. Using a 1:3 to 1:6 molar ratio of the vector DNA to the insert DNA.</li> <li>3. Do not try to amplify and sequence the clones when less than 100 transformed colonies obtained. Most of them contain false inserts.</li> <li>4. If the size of your insert DNA is about 1.6~1.8 kb, please perform colony PCR by specific primers to check the clone at first. If the clone is correct, then sequence it with M13 primers.</li> </ol>
High colony number, but high percentages of blue colonies with insert DNA.	<ol style="list-style-type: none"> <li>1. Improper ligation reaction.</li> <li>2. DNA is inserted, but it's not disrupting the expression of LacZ gene or the insertion generates in-frame fusion of the coding region to lacZ gene.</li> </ol>	<ol style="list-style-type: none"> <li>1. PCR DNA should be gel-purified or spin column-cleaned.</li> <li>2. Avoid multiple freeze-thaw and exposure to frequent temperature changes by preparing single-use aliquots of the ligase buffer.</li> <li>3. If a maximum of transformants is required, incubate the ligation reaction overnight at 4 °C.</li> <li>4. If multiple PCR products occur, target DNA should be gel-purified before ligation.</li> <li>5. Use colony PCR to screen both of the blue or white colonies.</li> </ol>
No blue colonies are present on the plate.	<ol style="list-style-type: none"> <li>1. Ampicillin is inactive.</li> <li>2. IPTG/X-Gal is insufficient or inactive.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check if ampicillin plates are prepared properly used within 1 month.</li> <li>2. Check if the IPTG/X-Gal are fresh and the IPTG/X-Gal plates should be stored in dark and cold place.</li> </ol>