


# NGS DNA Library Preparation

## With Single Stranded DNA

- **Simple Workflow**
  - Two reaction steps for library preparation
  - Simplified beads purification
- **Save Your Time**
  - Total time: 1~1.5 hours
  - Hands-on time: ~10 minutes
- **Application**
  - Single stranded DNA
  - Denatured DNA
  - Viral ssDNA



### Kit List

Product Name	Cat.#
NGS Single Stranded DNA Library Prep Kit (illumina platform, with non-index primers & PCR reagent)	30081
NGS Single Stranded DNA Library Prep Kit (illumina platform, with index primers & PCR reagent)	30082
NGS Single Stranded DNA Library Prep Kit (illumina platform, with unique dual index primers & PCR reagent)	30083