

IntelliXcap™

IntelliXcap is a next-generation capper and de-capper, engineered for increased throughput, ease of use and high system reliability, making the system ideal for any laboratory managing compound libraries or biological sample stores.

KEY FEATURES

Advanced Automation Means Increased Throughput

- ✔ IntelliXcap is extremely fast and works with multiple sample storage tube types from a range of tube manufacturers in 24, 48 and 96 format
- ✔ The only 96-format full rack capper and de-capper capable of de-capping a complete rack of 96 tubes in as little as 20 seconds
- ✔ Automation and high speed reduces sample handling time, increasing sample throughput workflow
- ✔ Ideal for medium to high-throughput laboratories managing compound libraries, biobanks or other biological samples

Compact, Modular Design Means Ultimate Flexibility


- ✔ IntelliXcap features the unique, fully automated, interchangeable IntelliXcartridge cap driver system
- ✔ IntelliXcartridge allows format change between different tube types in less than 2 minutes
- ✔ Simply by fitting the relevant IntelliXcartridge, IntelliXcap is compatible with all FluidX sample storage tubes, as well as sample storage tubes from Thermo Matrix, Thermo Nunc, LVL, Micronic
- ✔ IntelliXcartridge is available for both internal and external thread screw-cap tubes
- ✔ No need for specialist intervention to change cartridge drivers
- ✔ Control is via a quick installation, easy-to-use touch panel

Easy to Use, Easy to Integrate

- ✔ Quick installation, easy-to-use touch panels bring instant control
- ✔ Wide operating temperature range of 5°C to 40°C
- ✔ Electric torque control means less persistent wear on cap drivers
- ✔ Easily integrated, sample rack stage extension facilitates fully automated workflow integration

Weights and Dimensions

	Weight (kg)	Height (mm)	Width (mm)	Depth (mm)
IntelliXcap 96	22	320	256	468
IntelliXcap 48	28	386.2	256	464.4
IntelliXcap 24	27	386.2	256	464.4



Automated Storage Systems

Cryopreservation & Cold Chain Solutions

Informatics & Technical Solutions

Sample Storage, Lab Services & Transport

Sample Consumables & Instruments

Ordering Information

46-8012	FluidX IntelliXcap 96 96-Format Screw Cap Tube Rack Decapper/Capper with automatic glide rail for integration, includes 1 cartridge of choice, please quote tube type when ordering. Suitable for all FluidX tubes and almost all other brands with IntelliXcartridges
46-8011	FluidX IntelliXcap 48 48-Format Screw Top Tube Rack Decapper/Capper with automatic glide rail for integration, includes 1 cartridge of choice, please quote tube type when ordering. Suitable for all FluidX tubes and almost all other brands with IntelliXcartridges
46-8010	FluidX IntelliXcap 24 24-Format Screw Top Tube Rack Decapper/Capper with automatic glide rail for integration, includes 1 cartridge of choice, please quote tube type when ordering. Suitable for all FluidX tubes and almost all other brands with IntelliXcartridges

Ordering Information - IntelliXCap Cartridges: 96 format

48-8013-01	FluidX IntelliXCartridge 96-format Cap Driver Cartridge , with 96 individual cap drivers, 1 cartridge per case. Suitable for FluidX Internal Co-molded Thread
48-8013-02	FluidX IntelliXCartridge 96-format Cap Driver Cartridge , with 96 individual cap drivers, 1 cartridge per case. Suitable for FluidX External Thread
48-8013-03	FluidX IntelliXCartridge 96-format Cap Driver Cartridge , with 96 individual cap drivers, 1 cartridge per case. Suitable for FluidX Internal O-Ring Thread
48-8013-04	FluidX IntelliXCartridge 96-format Cap Driver Cartridge , with 96 individual cap drivers, 1 cartridge per case. Suitable for FluidX Acoustic Tube Thread
48-8013-05	FluidX IntelliXCartridge 96-format Cap Driver Cartridge , with 96 individual cap drivers, 1 cartridge per case. Suitable for Matrix 200ul Low Profile Internal Thread
48-8013-06	FluidX IntelliXCartridge 96-format Cap Driver Cartridge , with 96 individual cap drivers, 1 cartridge per case. Suitable for Thermo-Nunc Cryobank Internal Thread
48-8013-07	FluidX IntelliXCartridge 96-format Cap Driver Cartridge , with 96 individual cap drivers, 1 cartridge per case. Suitable for Micronic Internal Thread
48-8013-09	FluidX IntelliXCartridge 96-format Cap Driver Cartridge , with 96 individual cap drivers, 1 cartridge per case. Suitable for LVL Technologies Internal Thread
48-8013-11	FluidX IntelliXCartridge 96-format Cap Driver Cartridge , with 96 individual cap drivers, 1 cartridge per case. Suitable for Thermo/Matrix Internal Thread
48-8013-12	FluidX IntelliXCartridge 96-format Cap Driver Cartridge , with 96 individual cap drivers, 1 cartridge per case. Suitable for FluidX 0.2ml Tube External Thread

Ordering Information - IntelliXCap Cartridges: 48 format

48-8015-01	FluidX IntelliXCartridge 48-format Cap Driver Cartridge , with 48 individual cap drivers, 1 cartridge per case. Suitable for FluidX 1.5ml, 1.9ml, 3.8ml. Next-Gen Jacket and 2.0ml Jacket Automation-friendly External Thread
48-8015-03	FluidX IntelliXCartridge 48-format Cap Driver Cartridge , with 48 individual cap drivers, 1 cartridge per case. Suitable for FluidX 1.8ml Non-Jacket and Greiner 2.0 ml External Thread

Ordering Information - IntelliXCap Cartridges: 24 format

48-8017-01	FluidX IntelliXCartridge 24-format Cap Driver Cartridge , with 24 individual cap drivers, 1 cartridge per case. Suitable for FluidX 5ml Jacket and 7.6ml. Next-Gen Jacket Automation-friendly External Thread
48-8017-03	FluidX IntelliXCartridge 24-format Cap Driver Cartridge , with 24 individual cap drivers, 1 cartridge per case. Suitable for Micronic 6ml Internal Thread

Automated
Storage Systems

Cryopreservation &
Cold Chain Solutions

Informatics &
Technical Solutions

Sample Storage,
Lab Services & Transport

Sample Consumables
& Instruments



PROTEIGENE
Tél: 02.32.64.45.45 - Fax: 02.32.64.30.72
Email: info@proteigene.com - www.proteigene.com